

Information and criteria for buying Hyabell products

Hyabell filler products are 4th generation hyaluronic acid exclusively intended for use by professional aesthetic / medial technicians who master injection techniques for advanced dermal and subcutaneous fillers. Injecting of Hyabell requires specific knowledge and competence in the use, composition, indications, contraindications, incompatibilities and knowledge to treat side effects.

The new Hyabell generation of hyaluronic acid fillers have high viscoelastic properties. Elasticity and viscosity is directly proportional to the volume increasing and elevating properties Hyabell have in the skin tissue. Hyabell fillers have higher elastic modulus G and viscous G than any other cross-linked BDDE HA (hyaluronic acid) fillers. Injection results are immediately visible, and feel like a natural part of the skin tissue. The effect is long-term (from 6 up to 14 months).

Hyabell products are made of non-animal hyaluronic acid and contains 0.3% lidocaine. They are classified under medical equipment and designed to be injected by an experienced physician, especially in advanced injection techniques for filling in wrinkles and kontourering in the dermis and subcutaneous tissue.

In certain sensitive areas such as the lips, it is recommended to provide a local nerve blocking anesthesia. The doctor must master a number of established injection techniques for optimum results.

To order Hyabell products, we require a HPR number from the person placing the order. Responsible injectors must have certificate / diploma in advanced filler injections and documentation showing received practical training. Please contact us for an order form.

Scandinavian Skincare Systems UK
56 High Street, Broadway, Worcestershire WR12 7DT, England
+44 (0) 7711 086227

info@scandinavianskincareystems.com
www.scandinavianskincareystems.com